

Siena

Twin rail for curved stairs



Stannah

Make a step change

If a Stannah stairlift could make life easier, why wait?

Life is better when you have the energy to do the things you want to do, like enjoying a visit from the grandchildren or entertaining friends. Our experience has taught us that the sooner a stairlift is installed, the sooner you will be able to start enjoying the things you love doing. Just think how a Stannah stairlift could make your day-to-day life easier.

A stairlift that's right for you and your home.

We know that not all homes are the same and that curved staircases can come in different shapes and sizes. Our twin tube rail can be the ideal choice when you need a stairlift that will follow your staircase perfectly.



Inside rail

- Our stairlifts can accommodate different types of curved staircases.
- We can install our rails on either side of the staircase as shown in the photos below.
- Our slimline rails and fold-away chairs make the best use of available space and leave enough room for others to use the staircase.
- All our rails are tailor-made, so they fit your staircase perfectly.
- Contrary to popular belief our rails are fitted to the staircase and not the wall - ensuring minimal mess and bother.
- Installation is quick and easy.
- Our stairlifts can navigate all kinds of corners and landings with a range of rail solutions.



Outside rail

A company you can trust

We are a British, family-run and owned company and have been for five generations.

The current managing directors are the great, great grandchildren of Joseph Stannah, who started the company in the 1860s. His engineering principles and strong family values run as true today as they always have.

Now, over 150 years later, we have manufactured more than 700,000 stairlifts and installed products in over 40 countries.

Siena: Stylish, sleek and very comfortable

Our Siena stairlift combines fresh, contemporary styling with the best in stairlift engineering.

The simple contours along with five classic colour options make the Siena a smart and practical addition to any home.

The Siena comes with a twin rail that make the best use of space and with easy fold - just fold the arms then the seat base and footrest and the chair neatly folds away - so you can continue to use the stairs as usual if and when you want to.

The Siena also comes in two chair widths allowing you to pick the one most comfortable for you, your home and your life.



The Siena's slimline design allows others to still use the stairs

Siena key features

1. Footrest fold

To store the stairlift away, you can fold away the arms and seat and use the lever to raise the footrest. There's no awkward bending down, which is particularly useful at the top of the stairs.

Or for the ultimate ease-of-use opt for our powered footrest button.

2. Swivel seat

Using the levers enables you to swivel the seat easily, making it simpler to get off the chair. Alternatively, you can choose to have a powered swivel seat, which turns the chair to face the landing automatically.

3. Seatbelt

Choose a seatbelt to suit your style and buckle up for added security (retractable seatbelt shown – optional extra).



4. Safety edges

In the event that something obstructs the stairlift, sensors on the edge of the carriage and footrest will automatically bring the stairlift to a safe stop.

5. Easy-to-use controls

Effortlessly ride the Siena using our simple, intuitive controls that have been designed for those with less dexterity. They work with the lightest of touches, and can be operated using your hand, arm or wrist.

6. Wall controls

Ideal for when you need to summon your stairlift from another floor, or perhaps send it up or down to someone else who wants to use it.

Siena key benefits

Effortless and simple. This sums up beautifully just how much thought and care we have designed into each and every function of our Siena stairlift. From the safe and easy-to-use seatbelts, to stairlifts that fit on either side of the stairs and on all types of staircases, we have made your Siena stairlift simple, safe and reliable.



Seatbelt can be fastened with one hand



The Siena folds neatly away when not in use



Sensors detect any obstructions on your stairs

“We are delighted with the stairlift. My wife enjoys the ride. It gives her independence. She finds it very easy to use”

Mr A, Cirencester



The easy way to climb your stairs

Everything about a Stannah stairlift is made to make managing your stairs as easy, comfortable and safe as it can be, so your life changes for the better the moment you begin to use it.

From the simple controls at your fingertips to reliable operation and a flexible and compact design, we have done all the hard work so you don't have to. Just settle into the comfortable seat and sit back, as you travel up or down your stairs with ease.

We really have taken care of every detail to ensure your Siena stairlift is as effortless to use as it is easy to operate.

Step 1: Simply unfold the arms and seat base and then use the handy lever to lower the footrest.

Step 2: When the seat base and footrest are down, sit on your stairlift and make yourself comfortable in the

seat. Then pull across and fasten the easy reach seatbelt and you're now ready to go.

Step 3: A gentle movement of the directional control is all it takes for you to travel effortlessly and safely up or down your stairs.

Step 4: When you have arrived at the top of your stairs, you can then use your

swivel lever or optional powered swivel to rotate you onto your landing ready to get off your Siena in a safe position.

Step 5: When you are not using your stairlift, it can remain in position at the top of your stairs to act as a safety barrier, or, it can be folded away to allow space for other stair users.

“Now we have our Stannah stairlift, we no longer have to face an unwanted move into single storey accommodation”

Mr & Mrs H, Worcestershire



Tailored to you

Your comfort and safety mean everything to Stannah. This is why we have designed the Siena for complete peace-of-mind, giving it a range of helpful features that make it safe, easy and enjoyable to use.

Seatbelt choices



Lapbelt



Retractable seatbelt



Diagonal seatbelt



Five point harness

Control choices and control position

- Left hand
- Right hand



- Light-touch control



- Rudder style control

Folding footrest



- Powered footrest button

Wall controls



- Wall controls

Manual swivel



- Manual swivel seat

Powered swivel




- Powered swivel seat

Siena upholstery and rail colour choices


The Siena stairlift combines fresh, contemporary styling with the best in stairlift engineering. You can select one of our five classic faux leather collection or from one of ten colours of our quilted faux leather collection to blend in with your home decor.

Faux leather collection

A classic range of vinyl finishes to blend in with your home décor



Walnut



Ivory



Pine



Cranberry



Midnight

Quilted faux leather collection

An embossed stitched diamond faux leather range that comes in ten different colourways



Black



Maroon



Cocoa



Damson



Navy Blue




Green



Grey



Sand



Duck Egg Blue



Teal

Seat width options




Narrow seat



Standard seat

Rail colour options



Pearl (Standard)




Nut Brown



Jet Black



Window Grey



Beige



Ochre Brown

The pictures within our brochures are for guidance only. Although we have tried to display the colours accurately, we cannot guarantee that the colours in the pictures accurately reflect the colour of the actual stairlifts. Your stairlift may be slightly different from these pictures.

A professional installation

When you have chosen the Siena that suits you best, a Stannah trained engineer will visit your home and professionally install your new stairlift in a matter of hours, leaving your home just as they found it.





**Dolphin Stairlifts
(East Anglia Ltd)**
25 Estuary Drive,
Felixstowe, Suffolk, IP11 9TL

01394 277100
office@dolphinliftsanglia.com

Stannah

AUTHORISED DEALER

Our company policy is one of continuing product improvement. We therefore reserve the right to alter the specifications and to change and improve design at any time.

© Stannah Stairlifts Ltd 2020. Not to be reproduced wholly or in part without the written permission of Stannah Stairlifts Ltd. August 2021 DS